

可棄式注射過濾頭 Disposable Syringe Filter Units



· 小量氣體及液態過濾頭(DISMIC)

直徑 (mm)	孔徑 (μm)	盒入數	未滅菌	已滅菌	濾膜材質
3	0.20	100	03CP020AN	03CP020AS	醋酸纖維 Cellulose Acetate
	0.45	100	03CP045AN	03CP045AS	
	0.50	100	03JP050AN	-	疏水性 PTFE, Hydrophobic
13	0.20	100	13CP020AN	13CP020AS	醋酸纖維 Cellulose Acetate
	0.45	100	13CP045AN	13CP045AS	
	0.20	100	13HP020AN	-	親水性 PTFE, Hydrophilic
	0.45	100	13HP045AN	-	
	0.20	50	13HP020CN	-	親水性(離子層析) PTFE Hydrophilic for Ion Chromatography
	0.45	50	13HP045CN	-	
	0.20	100	13JP020AN	-	疏水性 PTFE, Hydrophobic
	0.50	100	13JP050AN	-	
	0.20	100	13NP020AN	-	尼龍 NYLON
25	0.20	50	25AS020AN	25AS020AS	混合纖維 MCE
	0.45	50	25AS045AN	25AS045AS	
	0.20	50	25CS020AN	25CS020AS	醋酸纖維 Cellulose Acetate
	0.45	50	25CS045AN	25CS045AS	
	0.80	50	25CS080AN	25CS080AS	
	0.20	100	25HP020AN	-	親水性 PTFE, Hydrophilic
	0.45	100	25HP045AN	-	
	0.20	50	25JP020AN	-	疏水性 PTFE, Hydrophobic
	0.50	50	25JP050AN	-	
	0.20	100	25NP020AN	-	尼龍 NYLON
	0.45	100	25NP045AN	-	
	0.20	50	-	25SS020RS	聚醚磺 PES
0.50	50	-	25SS045RS		

· 大量氣體及液態過濾頭(LABODISC)

直徑 (mm)	孔徑 (μm)	盒入數	未滅菌	已滅菌	濾膜材質
50	0.20	10	50CP020AN	50CP020AS	醋酸纖維 Cellulose Acetate
	0.45	10	50CP045AN	50CP045AS	
	0.20	10	50JP020AN	-	疏水性 PTFE, Hydrophobic
	0.50	10	50JP050AN	-	

◆ 產品編號

