

# Magnetic Stirrer

## MA100/300

Revolution(50/60Hz)	100~1,500rpm	400~1,550rpm	Stirring capacity	50~1,000ml	100~3,000ml
	MA100	MA300		MA100	MA300

### ■ MA100

- Small size
- Loaded with DC motor electronic control

### ■ MA300

- Long life
- Max. 3L stirring



MA100

MA300

### ■ Specifications

Model	MA100	MA300
Stirring stage size*	W77×D135mm *excludes liquid drop protection part	W178×D165mm *excludes liquid drop protection part
Stirring capacity	50 to 1,000ml	100 to 3,000ml
Revolution	Approx. 100 to 1,500rpm	Approx. 400 to 1,550rpm
Motor	DC motor, electronic control, 0.3W	AC motor, shaded pole type, 0.7W
Power (50/60Hz)	AC115V/AC220V Single phase with step-down transformer	AC115V/AC220V Single phase with step-down transformer
External dimensions	W84×D150×H60mm	W181×D197×H83mm
Weight	Approx. 1.1kg	Approx. 1.6kg
Accessory	Stirring Bar 30mm, 1pc.	

# Magnetic Stirrer

## MG120/600

Revolution(50/60Hz)	100~1,500rpm	200~1,500rpm	Stirring capacity	5~1,500ml×12 pcs.	50~2,000ml×6pcs.
	MG120	MG600		MG120	MG600

### ■ MG120

- Simultaneous Stirring up to 12 samples

### ■ MG600

- Individual Revolution Control for 6 points
- Chemical-resistant plate with ceramic coating



MG120

MG600

### ■ Specifications

Model	MG120	MG600
Stirrer	12 pcs. Interlock type	With individual revolution control for 6 points
Plate material	Ceramic coating	
Plate size	W295×D198 mm	W449×D300 mm
Stirring capacity	50 to 100ml×12 pcs. (Max.)	50 to 2,000ml×6pcs.
Revolution	Approx. 100 to 1,500rpm.	Approx. 200 to 1,500 rpm
Motor	DC motor×12pcs. Interlock electronic control	DC motor 0.7W×6pcs.
Power (50/60Hz)	AC115V/AC220V Single phase with step-down transformer	AC115V/AC220V Single phase with step-down transformer
External size	W303×D234×H55 mm	W456×D335×H90 mm
Weight	Approx. 4kg	Approx. 7.5kg
Accessory	Stirrer bar 25mm 12 pcs	Stirrer bar 30mm 6 pcs