

Loop Cinerator

SL-21

Internal sterilizing temp. 800~850°C

Time to reach sterilizing temp. 10min



■ Features

- Pt. Ear, Needle is sterilized in combustion tube.
- Safe operation through heater sterilization.
- Good heat-resistance.
- Loaded with quartz glass which heats evenly.
- Efficient for big volume and for continuous sterilizing.
-

■ Specifications

Model	SL-21
Internal sterilizing temp.	800~850°C
Time to reach sterilizing temp.	10min.
Cinerating port diameter	15mm
Heater	130W
Power	AC100V, 50/60Hz, 1.5A
External dimensions	W95 x D170 x H200mm
Weight	1.5kg

Flexible angles



Pt Ear can be placed in holder holes



Consumable

