

過濾器 Vacuum Filtration



KG-25



KG-47

- 微生物過濾，懸浮液中細菌，沈澱物等等微粒收集與分析

產品	KG-25	KGS-25
編號	17311200	17311300
Support Type	Sintered glass	Stainless steel screen

產品	KG-47	KGS-47
編號	17311400	17311500
Support Type	Sintered glass	Stainless steel screen



KG-90



KGS-90

產品	KG-90	KGS-90
編號	17312100	17312200
Support Type	Sintered glass	Stainless steel screen
備品更換		
1.Funnel, 1,100 mL	19312101	
2.Clamp	19312102	
3.Support screen (SUS316)	-	19312204
4.Base	19312103	19312203
5.#8A Stopper	19311404	




KP-47S



KSF-47

產品	KP-47S	產品	KSF-47
編號	43301010	編號	17312600
備品更換		備品更換	
1.Rubber cap	44501010	1.Funnel	-
2.Funnel cover	44501001	2.Support screen	19311504
3.Funnel	44501002	3.PTFE gasket	19311505
4.Funnel O-ring	44501008	4.Base	-
5.Support screen	44501004	5.#8A Stopper	19311404
6.Base	44501003		
7.#8B Stopper	19311008		

真空過濾器組 Vacuum Manifolds



KMP-3-6

- 抽氣式PVC過濾漏斗架，質輕、耐酸性
- 雙向閥由PTFE活塞與PVC主體組成

型式	編號	L x W x H cm	中心點間距 cm	重量 kg	
KMP-3	3 station	17313400	45.1 x 12 x 17.8	14	1.5
KMP-6	6 station	17313600	87.6 x 12 x 17.8	14	2.5



KM-3-6

- 抽氣式不鏽鋼過濾漏斗架，可高溫高壓滅菌
- 易沖洗，用做無菌測試

型式	編號	L x W x H cm	中心點間距 cm	重量 kg	
KM-3	3 station	17313000	45.7 x 12 x 17.8	14	5.0
KM-6	6 station	17313200	72.4 x 12 x 17.8	14	8.0